


Page
115-117

- Rubber Dam Clamps



Rubber Dam Clamp Forceps

www.surgiforce.com



RUBBER DAM FRAME

SF-5-5382
Young
105x95 Inox


SF-5-5383
Ivory
160 mm

SF-5-5384
Ainsworth
170 mm

SF-5-5385
Ainsworth
170 mm

SF-5-5386
Mershan
140 mm

CLAMPS (PLIERS)



SF-5-5387
Ivory light weight
170 mm

SF-5-5388
Brewer
175 mm

SF-5-5389
Brewer
175 mm

SF-5-5390
White
175 mm

SF-5-5391
Ivory
170 mm



SF-5-5392
Winged
Labial Clamp



SF-5-5393
Winged
Edged Labial Clamp



SF-5-5394
Winged
Cervical Clamp



SF-5-5395
Winged
Labial Clamp



SF-5-5396
Winged



SF-5-5397
Winged



SF-5-5398
Wingless

For Small Teeth, Elongated or With Narrow Collars



SF-5-5399
Winged



SF-5-5400
Winged



SF-5-5401
Wingless

For General Use With Adjustable Jaws According To Gingival Line Shapes



SF-5-5402
Winged



SF-5-5403
Edged, Winged



SF-5-5404
Wingless

For Premolars, With Elongated Collar



SF-5-5405
Winged



SF-5-5406
Edged, Winged



SF-5-5407
Wingless

For Premolars, With Elongated Collar



SF-5-5408
Wingless



SF-5-5409
Wingless

Rubber Dam Clamps

www.surgiforce.com



SF-5-5410
Winged



SF-5-5411
Wingless

For Small Lower Molars



SF-5-5412
Winged



SF-5-5413
Wingless

For Small Lower Molars With Round clamps



SF-5-5414
Winged



SF-5-5415
Wingless

For Square Lower Molars



SF-5-5416
Winged



SF-5-5417
Wingless

For Round Upper Molars



SF-5-5418
Winged



SF-5-5419
Wingless

For Molar Roots



SF-5-5420
Winged



SF-5-5421
Wingless

For Partially Broken Molars



SF-5-5422
Winged



SF-5-5423
Edged Winged



SF-5-5424
Wingless



SF-5-5425
Winged



SF-5-5426
Wingless

For Round Molars, Partially Broken, Undersized or Inclined



SF-5-5427
Wingless



SF-5-5428
Wingless



SF-5-5429
Winged



SF-5-5430
Winged



Surgiforce Instruments



Opposite Surgical Association,
Malkay Kalan, Marala Road,
Sialkot-51050 Pakistan.



+92 523 559131



info@surgiforce.com
sales@surgiforce.com



www.surgiforce.com



Our Certifications:

ISO 9001

CE Certified

EN 13485:2012

ISO 7153-1

ISO 7151

ISO 17664:2004

ISO 13402

ISO 15223

