



Surgiforce Instruments



DIAGNOSTIC instruments
Catalogue

www.surgiforce.com



Newest LED technology
One bulb is sufficient

Blade Size:

- No. 1
- No. 2
- No. 3
- No. 4



LED BULB

F.O. Blades:

- Mat Polished Stainless Steel
- Prevents Back Reflections
- Integrated Fiber Optics
- MEGALIGHT
- \varnothing 3-4 mm, 5500 Single Fibbers Minimum
- Can Be Treated In Autoclave
- 2.5, 3.5 V Battery Power
- Require 2 x C Size Batteries

FLAPLIGHT FIBER OPTIC LARYNGOSCOPE



Newest LED technology
One bulb is sufficient

Blade Size:

- Mac. 3
- Mac. 4



LED BULB

F.O. Blades:

- Mat Polished Stainless Steel
- Prevents Back Reflections
- Replacable Light Guide
- \varnothing 3-4 mm, 5500 Single Fibbers Minimum
- Can Be Treated In Autoclave
- 2.5, 3.5 V Battery Power
- Require 2 x C Size Batteries
- Swivelling blade tip facilitates difficult incubation.
- The Flaplight is a Macintosh laryngoscope blade with a tip which is adjustable through 70%



Newest LED technology
One bulb is sufficient

Blade Size:

- No. 00
- No. 0
- No. 1
- No. 2
- No. 3
- No. 4



LED BULB

F.O. Blades:

- Mat Polished Stainless Steel
- Prevents Back Reflections
- Integrated Fiber Optics
- MACINTOSH
- Ø3-4 mm, 5500 Single Fibbers Minimum
- Can Be Treated In Autoclave
- 2.5, 3.5 V Battery Power
- Require 2 x C Size Batteries

MILLER FIBER OPTIC LARYNGOSCOPE



Newest LED technology
One bulb is sufficient

Blade Size:

- No. 00
- No. 0
- No. 1
- No. 2
- No. 3
- No. 4



LED BULB

F.O. Blades:

- Mat Polished Stainless Steel
- Prevents Back Reflections
- Integrated Fiber Optics
- MILLER
- Ø3-4 mm, 5500 Single Fibbers Minimum
- Can Be Treated In Autoclave
- 2.5, 3.5 V Battery Power
- Require: 2 x C Size Batteries



LED BULB

Conventional Blades:

- Mat Polished Stainless Steel
- Prevents Back Reflections
- Can Be Sterilized
- Compatible According To DIN ISO 7376

Blade Size:

- No. 00
- No. 0
- No. 1
- No. 2
- No. 3
- No. 4

- Halogen and LED Bulb 2.5 V
- Require: 2 x C Size Batteries

MILLER CONVENTIONAL LARYNGOSCOPE



Blade Size:


- No. 00
- No. 0
- No. 1
- No. 2
- No. 3
- No. 4



LED BULB

Conventional Blades:

- Mat Polished Stainless Steel
- Prevents Back Reflections
- Can Be Sterilized
- Compatible According To DIN ISO 7376
- Halogen and LED Bulb 2.5 V
- Require: 2 x C Size Batteries



Blade Size:


- No. 1
- No. 2
- No. 3
- No. 4

Flange with an inter molded durable plastic locking mechanism.
The light source utilizes a crystal-clear plastic tubular fiber optic bundle style designed to reflect bright white light from the handle to the distal tip of the blade.

Economy Conventional Blades:

- Mat Polished Stainless Steel
- Prevents Back Reflections
- Can Be Sterilized
- Compatible According To DIN ISO 7376
- Halogen and LED Bulb 2.5 V
- Require: 2 x C Size Batteries

MILLER ECONOMY CONVENTIONAL LARYNGOSCOPE



Blade Size:

- No. 1
- No. 2
- No. 3
- No. 4

Flange with an inter molded durable plastic locking mechanism.
The light source utilizes a crystal-clear plastic tubular fiber optic bundle style designed to reflect bright white light from the handle to the distal tip of the blade.

Economy Miller Blades:

- Mat Polished Stainless Steel
- Prevents Back Reflections
- Can Be Sterilized
- Compatible According To DIN ISO 7376
- Halogen and LED Bulb 2.5 V
- Require: 2 x C Size Batteries



Attachment clip with integrated on/off switch. Secure. Switches off automatically when returned to the pocket.



Conventional
-Standard illumination
-With 2.5 V vacuum lamp

Fiber Optic
-High quality fiber optics
-High power 2.5 V xenon
-2.5 V LED lamp

-10 single use mini speculum (2.5 and 4mm) with storage bottle packing in black pouch
-Require: 2x AA batteries



FIBER OPTIC ADULT OTOSCOPES



-Bayonet fitting
-Illumination light 2.5v LED bulb
-Metal handle C
-Plastic Head CPP
-Rheostat dimmer
-05 x Ø 2.5 mm
-05 x Ø 4.0 mm
-Require: 2 x C batteries



Parker F.O. Otoscope set high quality professional Metal head, chrome plated battery handle Otoscope with extra bright light 2.5v LED bulb illumination provided thought fiber optic system. High quality big swiveling viewing window with 3x magnification lens out put for pneumatic membrane test.

-10 single use mini speculum (Ø 2.5 and 4mm) with storage bottle packing in black pouch.
-Require: 2 x C batteries



- Bayonet connecting system head fitting ideal for pneumatic tympanic membrane test, 2.5v LED white light
- 3 speculum (size Ø 2.5 3.5 and 4.5mm)
- 3x High quality magnifying lens
- Available: Polish and Matt finish
- Require: 2 x C batteries



- 2.5v standard illumination LED bulb. Chrome plated handle and head 4x magnifying glasses.
- Reusable speculum size Ø 2.5 3.5 and 4mm include. Connecting for pneumatic Otoscope. Battery handle with rheostat for Adjustment of illumination.
- 3 reusable ear specula.
- Available: Polish and Matt finish.
- Require: 2 x C batteries.

DERMATOSCOPES



- Supplies: Soft Black pouch packing with extra 0 - 10 scale glass plate
- Require: 2 x C size batteries.
- Available: Polish and Matt finish



- With approx. 10 times magnification
- focus adjustment – 6 up to +3.5 dpt
- halogen bulb that can be dimmed
- High capacity 2.5 V halogen lamp
- large contact disc (25 mm) with scale setting



Fiber Optic Battery Handles



Conventional Battery Handles

CONVENTIONAL LAYRNGOSOPE SETS



Macintosh Conventional
1 Handle + 3 Blades



Miller Conventional
1 Handle + 4 Blades



Miller Fiber Optic
2 Handles + 5 Blades



Macintosh Fiber Optic
1 Handle + 3 Blades

ENT DIAGNOSTIC SET



Set Boxes



ENT Set 2.5 V
With Accessories



Surgjforce Instruments



Opposite Surgical Association,
Malkay Kalan, Marala Road,
Sialkot-51050 Pakistan.



+92 523 559131



info@surgjforce.com
sales@surgjforce.com



www.surgjforce.com



Our Certifications:

ISO 9001

CE Certified

EN 13485:2012

ISO 7153-1

ISO 7151

ISO 17664:2004

ISO 13402

ISO 15223



No.3012922148