


Page
148–156

- Obstetrics
- Tracheotomy






SF-08-101
THOLE
for internal and external measurements, graduated in centimeters



SF-08-102
DE LEE
graduated in centimeters and inches

for internal measurements between symphysis and promontory



SF-08-103
COLLIN
graduated in centimeters




SF-08-104
MARTIN
graduated in centimeters




SF-08-105
COLLIN
graduated in centimeters

The additional straight arm can be replaced. The pelvimeter with straight arm is used for localisation of the portio vaginalis uteri and of vaginal metastasis.



SF-08-106
BREISKY
graduated in centimeters



SF-08-107
HERFF
puncturing forceps
225 mm




SF-08-108
IOWA
puncturing forceps
260 mm




1/3

SF-08-109
SIMPSON-BRAUN
305mm



1/3

SF-08-110
355 mm




1/3

SF-08-111
400 mm




1/3

SF-08-112
SIMPSON-BRAUN
305mm



1/3

SF-08-113
SIMPSON
235mm




1/3


SF-08-114
WRIGLEY
250 mm



BRAUN
Cranioclast
SF-08-115
420 mm




TARNIER
SF-08-116
400 mm




BRAUN
SF-08-117
310 mm



SMELLIE
SF-08-118
305 mm




BOER
SF-08-119
290 mm




SMELLIE
Perforator
SF-08-120
270 mm



**GREEN-ARMYTAGE
SF-08-121**
Caesarean section haemostatic
forceps, for edges of the transverse
incision in the lower segment
220 mm



**SANGER
SF-08-122**
275 mm



**SANGER
SF-08-123 - SF-08-124**
SF-08-123 270 mm
SF-08-124 290 mm




**GAUSS
SF-08-125**
Scalpel Flap
forceps
260 mm



**WILLETT
SF-08-126**
Scalpel Flap
forceps
185 mm




**SELLHEIM
SF-08-17**
320 mm




1/2


CUZZI
SF-08-128 - SF-08-130
 300 mm




SF-08-128
 17 mm



SF-08-129
 20 mm




SF-08-130
 23 mm



1/2

PESTALOZZA
SF-08-131
 305 mm




KANE
SF-08-132
 85 mm



COLLIN
SF-08-133
 85 mm






SF-08-141
KELLY
320mm




SF-08-142
SMELLIE
Perforator
270mm



SF-08-143
BRAUN
Cranioclast
420mm



SF-08-144
BOSE
165 mm



SF-08-145
BOSE
165 mm




SF-08-146 - SF-08-148
180 mm




SF-08-146
sharp




SF-08-147
blunt



SF-08-148
stumpf/blunt




SF-08-149
FINSEN
50 mm




SF-08-150
FINSEN
70 mm



SF-08-151
TROUSSEAU
120 mm



SF-08-152
LABORDE
125 mm



SF-08-153
LABORDE
strong pattern
135 mm



SF-08-154
Cleaning Brush
for tracheal tubes
285 mm



Surgiforce Instruments



Opposite Surgical Association,
Malkay Kalan, Marala Road,
Sialkot-51050 Pakistan.



+92 523 559131



info@surgiforce.com
sales@surgiforce.com



www.surgiforce.com



Our Certifications:

ISO 9001

CE Certified

EN 13485:2012

ISO 7153-1

ISO 7151

ISO 17664:2004

ISO 13402

ISO 15223



No.3012922148